OrthoPro® HyperEx Knee Orthosis

TREATMENT RATIONALE:

Provide knee hyperextension control for patients with genu recurvatum. A strategically placed air bladder is located on the hyperextension bar behind the knee. When inflated, the air bladder provides a gentle anterior directed rebound as the knee joint attempts to move into hyperextension.

FUNCTIONAL OBJECTIVES:

Provide hyperextension control allowing for greater independence in ADL'S. The HyperEx provides knee support and effective control of hyperextension with proper inflation of the air bladder.


ORTHOTIC TREATMENT AND FITTING:

1. Unfasten the thigh and shin straps on the lateral side of the brace.
2. With the patient seated and the leg fully extended, align the brace with the side bars centered on the leg on each side as follows:
 - a. The air nozzle and buckle attachments should be on the lateral side of the leg,
 - b. The air nozzle should be orientated at the top of the lateral side bar,
 - c. The hyperextension stop bar on the frame should be aligned behind the middle of the bend in the knee.
3. Tighten the thigh strap to comfort by securing the thigh strap onto itself over the thigh using the double-sided Velcro hook to hold the strap in place at the end of the strap. The end of the strap may be cut to size with scissors if too long. Do not cut too short to allow for potential weight gain.
4. Tighten the shin strap to comfort by securing the shin strap onto itself over the shin using the double-sided Velcro hook to hold the strap in place at the end of the strap. The end of the strap may be cut to size with scissors if too long. Do not cut too short to allow for potential weight gain.
5. Have the patient stand. Adjust the thigh and shin straps to comfort.
6. Adjust the hyperextension strap as needed to optimally stop unwanted hyperextension or Recurvatum of the knee joint. The adjustable hyperextension strap is adjusted by tightening or loosening the straps on either side of the brace at mid knee and securing them back on the adhesive hook located on the rigid hyperextension stop bar.
7. Using the air actuator, insert clear tip into the bladder valve to comfortably cushion the space between the hyperextension stop strap and the rigid extension stop frame bar.
8. Have the patient walk in the brace. Adjust the strapping and air bladder to ensure comfort and the desired Recurvatum stop action of the brace.
9. To remove air from the bladder, insert the black tip of the air actuator into the bladder valve.
10. Always remove air before removing the brace and re-inflate when applying the brace
11. Follow manufacturer's instructions for care of the device. Always inspect the device between applications to ensure the soft goods are properly in place, the device has not been badly soiled or would provide any other risk to the patient.

LAUNDRY INSTRUCTIONS

1. Brace and strap pads can be wiped clean with a damp cloth and/or antibacterial wipe.
2. If excessively soiled, remove pads from straps and brace (secured by hook and loop) and hand wash with mild detergent. Towel dry and lay flat until completely dry. **DO NOT PUT IN WASHER AND DRYER.** Reapply to brace when completely dry.

WARNING: The product should be fit by trained personnel. The product is designed for single patient use only in order to avoid cross contamination. Any substitution or removal of the product's parts voids the manufacturer's warranty. OCSI, Inc. will assume no liability if the above instructions are not followed.

Authorized Representative
MDSS GmbH, Schiffgraben 41
30175 Hanover, Germany

Ongoing Care Solutions, Inc
11721 US HWY 19 N
Clearwater, FL 33764
800-375-0207 (P) 727-525-1424(F)
www.ongoingcare.com

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Authorized Representative
MDSS GmbH, Shiffgraben 41
30175 Hanover, Germany

Ongoing Care Solutions, Inc
11721 US HWY 19 N
Clearwater, FL 33764
800-375-0207 (P) 727-525-1424(F)
www.ongoingcare.com