AirPro® Grip WHFO

TREATMENT RATIONALE:

To treat joint stiffness, or contractures of the wrist / hand / and fingers (including thumb); ideal for treating a fist hand. An inflatable cone shaped air bladder located on the palmar bar is carefully inserted deflated into a fist hand. Adjust the wrist, forearm support to the desired position. Using an inflation/deflation bulb, air is inserted into the air bladder to provide a gentle stretch on the affected joints of the fingers and thumb to reverse significant adaptive tissue shortening of the fingers and thumb. The wrist and forearm are supported at the desired angle by hand molding the long opponet of the brace to the desired angle. The orthotic treatment should be continued until function is restored to the affected hand.

FUNCTIONAL OBJECTIVES:

To effectively treat non-fixed contractures of the wrist, hand, and fingers. Ideal for therapeutic stretching of the fingers and thumb while supporting the wrist and forearm.


ORTHOTIC TREATMENT AND FITTING:

1. Use PROM or NeuroStretch™ submaximal passive stretching to point of noticeable resistance only (no discomfort) to passively stretch the affected joint capsule(s), connective tissue, tendons, and muscles. Concentrate on the wrist and thumb NeuroStretch™ locations prior to placing any extension force on the fingers or thumb.
2. Gently insert the deflated palmar bar into the palm of the hand. Make sure that the cone is centered as well as possible with each end of the cone extending from each side of the hand.
3. Secure the forearm strap first. The strap should be snug, but not too tight. Next, apply the wrist strap to comfort. When both the forearm strap and wrist strap have been secured, apply the finger strap upwards over the inflation nozzle and over the fingers. The slit in the finger strap should go over the tip of the small end of the palmar cone. Secure the finger strap back onto itself to maintain the position of the fingers in the device.
4. The inflation/deflation valve is on the larger end of the cone. Using the inflation/deflation bulb, place the clear end of the bulb into the valve and slowly inflate the palmar cone. As the cone is inflated, it will become more difficult to continue to squeeze the inflation bulb as air fills the air bladder. At this point, air pressure is providing a stretch on the affected fingers. Use caution at this point to NOT over inflate the air bladder causing unwanted pressure or over stretching of the fingers and thumb. The patient should experience a gentle stretching with no significant discomfort or pain. The black tip of the bulb is the deflation end of the inflation/deflation bulb, and can be used to reduce air pressure as needed.
5. The orthotic wrist base (long opponet) may be hand molded as needed. Provide approximately 5 to 10° of additional wrist stretch to maximize patient outcomes in treating wrist flexion.
6. After initial device application the affected tendon(s) should feel stretched with no indication of pain or discomfort. After 15 minutes of wear, softening or relaxation of the same tendon(s) indicates that the joint is predisposed to long effects therapeutic stretch.
7. Gradually increase wearing time to three hours or more to achieve Total End Range Time (TERT) that will provide long effects stretch. For a severely contracted hand, careful monitoring and gradual increase in wearing time is recommended.
8. Always deflate the palmar cone using the black tip of the inflation/deflation bulb prior to removing the device.
9. Determine wearing schedule based upon patient tolerance, therapy evaluation and physician's order.
10. Check for skin redness, pressure, and potential patient discomfort every two to three hours. Evaluate any red areas using the Blanch Test. Remove the orthotic device immediately if significant redness, pressure, or pain and discomfort are evident. ***If there is significant redness, pressure, or pain associated with device use, remove the device immediately.*** Discontinue use until the skin integrity or comfort issues are resolved. The device may require modification or the wearing schedule may be altered to eliminate potential skin integrity and comfort problems.
11. A significant wearing schedule disruption often requires a re-adaptation period. A gradual re-introduction and increased wearing endurance may be necessary. Note wearing schedule changes in the Patient's care plan.
12. Follow manufacturer's instructions for care of the orthotic device. Always inspect the device between applications to ensure the soft goods are properly in place, the device settings have not been altered, and the device has not been soiled or would provide any other risk to the patient prior to application.
13. Whenever the patient is not experiencing a gentle stretch sensation post application (or at least once a month), check the device settings for continued application of the desired amount of progressive extension therapeutic stretch.

Laundry Instructions:

1. Always remove soft cover from frame before washing. Remove the air bladder from the palmar cone.
2. Close hook and loop attachments and place in enclosed laundry bag.
3. Hand or machine wash, gentle cycle with mild detergent. No bleach or fabric softener.
4. No bleach or fabric softener.
5. Air dry.

WARNING: The product should be fit by trained personnel. The product is designed for single patient use only in order to avoid cross contamination. Any substitution or removal of the product's parts voids the manufacturer's warranty. OCSI, Inc. will assume no liability if the above instructions are not followed.

Authorized Representative
MDSS GmbH, Schiffgraben 41
30175 Hanover, Germany

Ongoing Care Solutions, Inc
11721 US HWY 19 N
Clearwater, FL 33764
800-375-0207 (P) 727-525-1424(F)
www.ongoingcare.com

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