DynaPro® Torticollis Orthosis

TREATMENT RATIONALE:

Treatment of lateral flexion of the cervical spine / Torticollis. The cervical spinal lateral flexion condition must be functional in nature / correctable (not structural) to be considered medically necessary.

FUNCTIONAL OBJECTIVES:

Improve Range of Motion and functional alignment of the head and neck. Reduce cervical flexion due to abnormal tone and spasticity. Support weak neck muscles to assist in functional alignment of the head and neck. Improve patient ability to swallow and interact with the social environment with improved head positioning.

ORTHOTIC TREATMENT AND FITTING:

1. Determine the correct size of DynaPro® Torticollis Orthosis. The device is ordered as Left or Right depending on the direction of lateral flexion of the head / neck to be treated.
2. Determine comfortable range of motion, resistance levels, and need for spinal immobilization or support through a cervical spinal assessment.
3. Modify the thermoplastic insert as needed to support the neck, bringing the cervical spine towards a more functional alignment. The head / neck portion of the device should be set approximately 15° of additional beyond significant resistance to stretch towards functional alignment to reverse Torticollis.
4. Apply the DynaPro Torticollis Orthosis on the patient's shoulder by bending the shoulder pad to conform to the patient's shoulder with the head / neck piece aligned properly at the mid line of the head and neck. Check for the thermoplastic head / neck upright to ensure the absence of any pressure on the ear or TMJ joint.
5. With the orthosis correctly aligned, bring the securing strap behind the back, under the opposing axilla, and over the chest. The end of the securing strap should be placed through the securing loop on the anterior side of the shoulder pad and attached back onto itself. The axilla pad should be centered under the axilla to provide padding from the strapping system.
6. After the strapping has been secured, check the device alignment and fit. Adjust the strapping if necessary.
7. Palpate the muscles on the shortened side of the neck immediately after device application. The muscles should be taut indicating that a stretch is being applied on the shortened muscles by the device.
8. Palpate the shortened muscles approximately 5 minutes after device application. There should be significant softening of the stretched muscles demonstrating that the patient has relaxed into the device. Adjust the device to lessen the stretch if necessary (no softening).
9. Allow the patient to stand or be seated and make any adjustments as necessary. Device should be worn while standing or sitting.
10. Determine wearing schedule per therapy evaluation and physician's order.
11. Always inspect the skin after removal to ensure the device applies no unwanted pressure.
12. Modify incrementally over time as necessary to incrementally improve cervical alignment to neutral

LAUNDRY INSTRUCTIONS:

1. Always remove soft cover from the orthotic frame before laundering.
2. Close hook and loop attachments and place in enclosed laundry bag.
3. Hand or machine wash, gentle cycle with mild detergent. No bleach or fabric softener.
4. Air dry.

INSERT:

1. The aluminum / foam shoulder pad is medical grade closed cell foam. An antimicrobial wipe can be used to wipe down the entire surface of the shoulder pad. Allow to air dry.
2. The thermoplastic head / cervical part of the device can be wiped with an antimicrobial wipe. The device should be allowed to air dry before placing the insert back into the soft cover.

WARNING: The product should be fit by trained personnel. The product is designed for single patient use only in order to avoid cross contamination. Any substitution or removal of the product's parts voids the manufacturer's warranty. OCSI Inc. will assume no liability if the above instructions are not followed.

MD

Authorized Representative
MDSS GmbH, Schiffgraben 41,
30175 Hanover, Germany

Ongoing Care Solutions, Inc
11721 US HWY 19 N
Clearwater, FL 33764
800-375-0207 (P) 727-525-1424(F)
www.ongoingcare.com

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