OrthoPro® ROM Knee Orthosis

TREATMENT RATIONALE:

Facilitate maintenance or rehabilitation of the knee in treatment of joint stiffness or non-fixed adaptive tissue shortening contractures of the knee joint. The orthosis ROM hinge can be set in 15° steps to control both flexion and extension motion in the brace. Ideal for providing protected ROM exercise in the brace post surgically or to limit motion to prevent further injury of the joint during rehabilitation. The OrthoPro™ ROM Knee can also be used at a static setting to protect and support the joint or to provide progressive extension orthotic therapy.

FUNCTIONAL OBJECTIVES:

Maintain or increase elbow Range of Motion, allowing for increased functional use of the affected joint. Functional use of the knee can be significantly improved by increasing ROM, allowing for greater independence in ADL'S. If function is not restored or significantly limited, the device should be used to maintain ROM as needed. Provide support of the knee joint during weight bearing and initial assisted ambulation.

ORTHOTIC TREATMENT AND FITTING:

1. Use PROM to passively stretch the affected joint capsule(s), connective tissue, tendons, and muscles.
2. Slowly and gently stretch the knee to the point of noticeable resistance only (no discomfort). Hold for a minute to allow the extension release of the affected joint.
3. Set the elbow orthotic hinges to provide the desired orthotic therapy. Follow the instructions for the ROM hinge setting. For a progressive extension orthotic therapy setting, add approximately 15° of extension to the comfortable end range for the flexion stop hinge setting. When the device is applied, a gentle stretch will be applied to the affected joint. Open the upper and bottom cuff flaps. The device label is always placed on the top cuff for ease in device orientation.
4. Holding the knee at its furthest point of extension, place the orthotic on the leg with the hinge at the midline of the knee joint.
5. Secure the top cuff flap over the top of the thigh.
6. Secure the bottom cuff flap over the top of the shin.
7. Place the kneecap directly over the patella.
8. Secure hook end of the kneecap straps across the kneecap in a "X" pattern to evenly distribute any pressure across the entire knee cap control pad.
9. Check all cuff flaps and the knee control pad for pressure. Two fingers should be able to be inserted between the flap and the skin. Loosen flaps if necessary.
10. Determine wearing schedule per therapy evaluation and physician's order.
11. Incrementally increase wearing time per patient tolerance and patient care plan up to a maximum of six hours on per shift. A minimum of three to four hours of wear daily after the adaptation period is recommended for best results.
12. Release and check for skin redness or pressure or patient discomfort every two to three hours at a minimum. Remove the orthotic device immediately if significant redness or pressure is evident. Notify the appropriate staff member(s) immediately and document any significant redness or signs of adverse pressure or shear. Discontinue device use until the skin integrity issues are resolved, and the device is modified, or the wearing schedule is altered to eliminate potential skin integrity problems.
13. Follow manufacturer's instructions for care of the orthotic device. Always inspect the device between applications to ensure the soft goods are properly in place, the device settings have not been altered, and the device has not been soiled or would not provide any other risk to the patient prior to application
14. Check device settings for continued application of the desired amount of extension at least once a month. Re-adjust the hinges to maintain 5° to 10° of extension beyond the point of resistance to stretch as needed for progressive extension orthotic therapy.


LAUNDRY INSTRUCTIONS

- Remove uprights.
- Fasten all hook and loop attachments on soft cover.
- Hand or machine wash, gentle cycle with mild detergent. DO NOT USE COMMERCIAL WASHERS OR HOT WATER.
- No bleach or fabric softener.
- Air dry.

HINGES:

- Clean with mild detergent
- Rinse with clean water
- Dry with soft cloth

WARNING: The product should be fit by trained personnel. The product is designed for single patient use only in order to avoid cross contamination. Any substitution or removal of the product's parts voids the manufacturer's warranty. OCSI, Inc. will assume no liability if the above instructions are not followed.

Authorized Representative
MDSS GmbH, Shiffgraben 41
30175 Hanover, Germany

 **Ongoing Care Solutions, Inc**
11721 US HWY 19 N
Clearwater, FL 33764
800-375-0207 (P) 727-525-1424(F)
www.ongoingcare.com

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Authorized Representative
MDSS Gmbh, Shiffgraben 41
30175 Hanover, Germany

Ongoing Care Solutions, Inc
11721 US HWY 19 N
Clearwater, FL 33764
800-375-0207 (P) 727-525-1424(F)
www.ongoingcare.com