SoftPro® II In-Bed Ankle / Foot Orthosis

TREATMENT RATIONALE:

Orthotic therapy for mild to moderate joint stiffness, contracture, or abnormal tone and spasticity of the ankle and foot. The orthotic device treats plantar flexion and inversion /eversion of the ankle / foot, and rotation of the hip. A removable adjustable strap is included with the device that can be used to treat inversion/eversion or be used as a dorsi-assist. Although a recumbent orthosis, the SoftPro® In-Bed II AFO has a non-skid pad sewn into the bottom to facilitate transfers without having to remove the AFO. The AFO is made of a flexible material, and provides a gentle dorsi-flexion stretch of the ankle / foot.

FUNCTIONAL OBJECTIVES:

Increase available range of motion, active and/or passive, at the ankle joint to:

- Improve / maintain functional alignment of the ankle / foot.
- Increase range of motion and mobility.
- Decrease risk of skin breakdown.
- Reduce / eliminate contracture related pain and discomfort.
- Prevent permanent deformity.

ORTHOTICE TREATMENT AND FITTING:

1. If necessary, heat mold the semi-rigid insert if the plantar flexion or inversion /eversion is greater than 15° from neutral. Mold to within 10° of comfortable end range stretch for plantar flexion, and at end range for inversion/eversion and to accommodate any rotation of the foot to significantly reduce any possibility of pressure on the ankle/foot with device use.
2. Open the closures so that the device is ready to apply after passive stretching of the ankle/foot.
3. With the patient in a supine position, passively stretch the ankle/foot. Flex the knee if necessary to 45° or as close to 45° as comfort will allow obtaining optimal positioning of the ankle/foot post stretch. Avoid a stretch reflex while passively stretching the joint. Passive stretch must address the plantar flexion, inversion /eversion, and rotation of the foot. Slowly and gently, use sub-maximal passive stretching to point of noticeable resistance only (no discomfort). Hold for a minute to allow the extension release of the affected joint(s).
4. Place the patient's foot in the AFO insuring that the heel of the foot is at the apex of the heel of the AFO.
5. Holding the AFO in place, bring the loop side of the foot portion of the liner over the top of the foot. Bring the opposite side of the foot liner over the top of the foot and secure the hook and loop closure. A finger should comfortably slip under both edges of the closure to ensure that the closure is not secured too tightly.
6. Bring the loop side of the calf portion of the liner over the calf. Bring the opposite side of the liner over the calf and secure the hook and loop closure. A finger should comfortably slip under both edges of the closure to ensure that the closure is not secured too tightly.
7. Place the anterior ankle pad over the top of the ankle and bring the unsecured end of the ankle strap through the open slot at the bottom of the heel, and secure the strap back onto itself. Adjust to a comfortable fit, secure enough to hold the ankle in place. The strap should be snug, but not tight. If needed, adjust both ends of the anterior strap if the strap needs to be centered to maintain a good closure.
8. Bring the hook closure at the top of foot liner portion of the liner over the anterior strap and attach the closure over the strap.
9. Check to insure the heel is floating free in the brace with clearance off of the plastic heel of the brace.
10. Move the hip control bar if needed to the left or right side to control or prevent internal or external rotation of the hip when the resident is in bed.
11. Attach the removable pull strap by placing the center of the strap under the top of the foot plate and securing the hook and loop. Take each end of the strap up through the D ring on the top back of the AFO, and secure the strap onto itself at the desired therapeutic tension. Tighten both straps to treat plantar flexion.
12. If using the pull strap, tension should initiate a GENTLE stretch of the ankle/foot in the corrective direction. There should be no pain or discomfort. The patient should experience a gentle stretching for several minutes. The stretching sensation should diminish and become more comfortable after several minutes of wear. White out or an indelible marker can be used to mark the pull straps on each side to delineate pull strap settings for caregivers.
13. Determine wearing schedule per therapy evaluation and physician's order.
14. Use the device in the recumbent position only. ***This device is not an ambulating AFO.***

IMPORTANT:

- **Do not apply the device if there is significant skin redness that would come in contact with the device**
- **Upon removal of device, closely inspect skin integrity. If any redness is present, discontinue device use until skin integrity issues are resolved and the orthosis is adjusted and/or the wearing schedules adjusted accordingly.**


CONTRA – INDICATIONS:

- An ankle/foot with grade three plus edema
- An ankylosed ankle or and ankle/foot that is broken or dislocated


LAUNDRY INSTRUCTIONS:

1. Always remove soft cover from frame before washing. Wipe frame with damp cloth or antibacterial wipe.
2. Close all hook and loop attachments on soft cover.
3. Hand or machine wash, gentle
4. Air or tumble dry low heat.
5. No bleach or fabric softener.
6. DO NOT USE COMERCIAL WASHER OR DRYER!

WARNING: The product should be fit by trained personnel. The product is designed for single patient use only in order to avoid cross contamination. Any substitution or removal of the product's parts voids the manufacturer's warranty. OCSI Inc. will assume no liability if the above instructions are not followed.

Authorized Representative
MDSS Gmbh, Schiffgraben 41
30175 Hanover, Germany

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
MD

Authorized Representative
MDSS GmbH, Schiffgraben 41
30175 Hanover, Germany

Ongoing Care Solutions, Inc
11721 US HWY 19 N
Clearwater, FL 33764
800-375-0207 (P) 727-525-1424(F)
www.ongoingcare.com