



GUARDIAN • BRACE •

REHABILITATOR™ KNEE MEASUREMENTS

OA REHABILITATOR™ SIZING CHART


SIZE	Circumference 6" above mid patella (B)	Circumference 4" below mid patella (D)
XSMALL	13" (33cm) – 16" (40.6cm)	9" (22.9cm) - 12" (30.5cm)
SMALL	16" (40.6cm) – 19" (48.3cm)	12" (30.5cm) - 15" (38.1cm)
MEDIUM	19" (48.3cm) – 22" (55.9cm)	15" (38.1cm) - 17" (43.2cm)
LARGE	22" (55.9cm) - 25" (63.5cm)	17" (43.2cm)- 20" (50.8cm)
XLARGE	25" (63.5cm) – 28" (71.1cm)	20" (50.8cm) - 23" (58.4cm)
XXLARGE	28" (71.1cm) - 32" (81.3cm)	23" (58.4cm) - 26" (66cm)

SPORT REHABILITATOR™ SIZING CHART

SIZE	Circumference 4 1/2" above mid patella (C)	Circumference 4" below mid patella (D)
XSMALL	13" (33cm) – 16" (40.6cm)	9" (22.9cm) - 12" (30.5cm)
SMALL	16" (40.6cm) – 19" (48.3cm)	12" (30.5cm) - 15" (38.1cm)
MEDIUM	19" (48.3cm) – 22" (55.9cm)	15" (38.1cm) - 17" (43.2cm)
LARGE	22" (55.9cm) - 25" (63.5cm)	17" (43.2cm)- 20" (50.8cm)
XLARGE	25" (63.5cm) – 28" (71.1cm)	20" (50.8cm) - 23" (58.4cm)
XXLARGE	28" (71.1cm) - 32" (81.3cm)	23" (58.4cm) - 26" (66cm)

SPORT EXT REHABILITATOR™ SIZING CHART

SIZE	Circumference 9" above mid patella (A)	Circumference 6" below mid patella (E)
XSMALL	15" (38.1cm) - 18" (45.7cm)	9" (22.9cm) - 12" (30.5cm)
SMALL	18" (45.7cm) - 21" (53.3cm)	12" (30.5cm) - 15" (38.1cm)
MEDIUM	21" (53.3cm) - 24" (61cm)	15" (38.1cm) - 17" (43.2cm)
LARGE	24" (61cm)- 28" (71.1cm)	17" (43.2cm)- 20" (50.8cm)
XLARGE	28" (71.1cm)- 31" (78.7cm)	20" (50.8cm) - 23" (58.4cm)



9" above mid Patella (A)

6" above mid Patella (B)

4 1/2" above mid Patella (C)

4" below mid Patella (D)

6" below mid Patella (E)

6551 43rd Street North, Unit 1402
Pinellas Park, FL 33781
(800) 375-0207